

LED PENLIGHT



LUMINDEX 3

The light has been manufactured with a super-high luminosity LED. The intensity is from 15,000 to 30,000 Lux at 5 mm. These characteristics enable the light to support diagnosis during a therapeutic or surgical intervention.

- Super-high luminosity LED with assists in inspecting caries, calculus, fractures, and root canal orifices.
- Can reduce the need and expense of diagnostic X-rays.
- Immediate feedback saves time and improves patient care.
- Simple HI-LO-OFF one-button, battery-powered operation.

Visual imaging fracture, caries and root canals



Contents

- LUMINDEX™ 3 handpiece x1
- 3 mm light guide rod x1
- AAA alkaline battery x2
- Disposable sleeve x20
- Instructions for use x1



Specification

Battery type	2 x AAA alkaline batteries
Color temperature	White LED, 6500 K
Radiant intensity	Low mode : 15,000 Lux above at 5 mm High mode : 30,000 Lux above at 5 mm
Dimensions	16 mmØ (max)x 135mm
Weight	59 g (including batteries)

Optical Accessories: Endo-Lux

The more Endo-Lux performs, the more posterior caries, cracks and fractures will be identified.



2402-0018 Endo fiber 1mm,40pcs/set



2402-0017 Endo-Lux
5.5 Ø 1 with 10 pcs Endo fiber